




Cytomass Monitor

Measure the mass of living
cells in real time





High-quality frequency and phase data are essential for precise mass measurements. This is why the Cytomass Monitor uses photo-excitation of its cantilever. Observe the excellent correlation of experimental data (black squares) and fitted curves (red lines).




Optical microscope integration and full temperature, gas, and humidity control allow continuous monitoring of cells over a long period of time (up to days) with optical techniques (e.g. DIC, fluorescence) in addition to mass measurements. Images in graph show a time series of Histone H2B-GFP-labelled HeLa cells attached and growing on the Cytomass Monitor's cantilever. Data courtesy Müller lab ETHZ.

Measuring and monitoring the mass of a single cell in physiological conditions

Cells tightly regulate their mass and volume during development and homeostasis. Though this regulation is fundamental to proper functioning of the cell, it is so far poorly understood, mainly because current tools can only reasonably determine the volume of spherical suspended cells (~1% population accuracy, ~20% single cell accuracy), not their mass and not under normal adherent growth conditions.

Together with Prof. Daniel J. Müller and Dr. David Martínez-Martín of the ETH Zurich, Nanosurf is developing the Cytomass Monitor. This cantilever-based instrument is able to measure the mass of even single adherent cells with 0.1% mass and 10-ms temporal resolution.



$$f_0 = \frac{1}{2\pi} \sqrt{\frac{k}{m^*}}$$

f_0 : Resonance freq.
 k : Spring constant
 m^* : Effective mass

Working principle

Key Features

- Ability to accurately and continuously measure (monitor) the mass of cells in physiological conditions (cell culture) over days with millisecond resolution
- Ability to measure the mass of adherent cells while being anchored to physiologically relevant substrates (e.g. collagen, fibronectin, vitronectin, laminin)
- Ability to measure the mass of a group of cells or of a small piece of tissue, allowing to study the influence of cell-cell interactions on mass changes
- Ability to measure the mass of suspended cells by immobilizing them
- Fully compatible with state-of-the-art cell biology optical microscopy techniques

Key applications

- Cell mass or volume regulation
- Cell migration
- Cell nutrition
- Cell division and cell cycle progression
- Fat cell storage and metabolism
- Viral infection-related processes
- Ion channel properties
- Drug screening for drugs targeting pathways linked to cell growth
- New therapies for cancer, aging, obesity, type 2 diabetes, neurodegeneration, and other diseases linked to a deregulation of growth control

Groundbreaking research

Interested in becoming a Cytomass Monitor pioneer?

Then stay ahead of your peers and contact us now at info@nanosurf.com

Nanosurf AG

Gräubernstrasse 12
4410 Liestal
Switzerland
+41 61 927 47 47 (phone)
+41 61 927 47 00 (fax)

Nanosurf GmbH

Rheinstrasse 5
63225 Langen
Germany
+49 6103 202 7163 (phone)
+49 6103 202 7182 (fax)

Nanosurf Inc.

300 Trade Center, Suite 5450
Woburn, MA 01801
United States of America
781 549 7361 (phone)
781 549 7366 (fax)

Nanosurf 中国

Nanosurf China, Shanghai
上海市天宝路578号 (200086)
飘鹰世纪大厦703室, 中国
+86 18621896399 (电话)
+86 21 5512 7698 (传真)