

Elasticity measurements with colloidal probes

Advantages of FluidFM compared to classical AFM methods

AFM can be used to determine the mechanical properties of cells, tissue, and substrates. These mechanical properties have shown to be important for fundamental cellular processes like migration^[1] and differentiation^[2], as well as in the characterization and understanding of cancer^[3].

Although colloidal probes are commercially available, it is still common practice to glue beads on the cantilever. This is a precision job and can only be done one bead at a time. Coating has to be applied with the glued bead on the cantilever. With FluidFM, gluing of beads can be circumvented^[4]: pre-coated beads are mounted to the hollow cantilever shortly before the experiment by application of underpressure. Upon contamination, a bead can be immediately exchanged and measurements can continue with the same cantilever.





Conventional AFM methods

Preparation:

- Time-consuming and complex gluing routine needs preparation and experience
- Coating procedure on complete cantilever, cumbersome to parallelize

Experiment:

- A simple replacement of the bead after contamination or for experiments using differently coated beads is impossible



Colloidal spectroscopy of HeLa cells using FluidFM. (A) SEM image of cantilever. (B) Optical view of cantilever, colloidal beads, and cells in the Cytosurge control software. (C) Map of Young's modulus (scan size: 60 μm; blue area is glass). (D) Approach curves (force in nN versus tip-sample distance in μm) for the locations indicated with a star in C: glass (blue), harder cell areas (orange) and softer cell areas (black).

Advantages of FluidFM

- Handling — Pre-coated beads and easy mounting
- Experiment design — Same cantilever and different beads
- Time — Fast mounting and bead exchange
- Easy — Intuitive touchscreen-based workflow environment

[1] Park et al. 2005, *Cell Motility and Local Viscoelasticity of Fibroblasts*, Biophys. J. 89, 4330–4342

[2] Engler et al. 2006, *Matrix Elasticity Directs Stem Cell Lineage Specification*, Cell 126, 677–689


[3] Plodinec et al. 2012, *The nanomechanical signature of breast cancer*, Nat. Nanotechnol. 7, 757–765

[4] Dörig et al. 2013, *Exchangable colloidal AFM probes for the quantification of irreversible & long-term interactions*, Biophys. J. 105, 463–472

Typical elasticity measurement series with self-attached micro-spheres*.

FluidFM and conventional AFM compared.

Preparation: Conventional AFM method





FluidFM



Batch processing of multiple beads: a large reservoir of equally coated beads is accessible.

Experiment:



* The above schematic represents a simplification of the steps actually performed for each method. They may vary with the goals for each experiment and method. Although not ideal scientifically, conventional AFM may be conducted differently to overcome some of the inherent limitations of the method, and FluidFM could be used in a more comparable "classical" way. Nevertheless, changing the bead without changing the cantilever remains inaccessible to conventional AFM, and FluidFM will always provide more experimental freedom and faster performance.